OMRON

# GS CUFF GS カフ

SS (HXA-GCUFF-SSLA / HXA-GCUFF-SSLB)  
S (HXA-GCUFF-SLA / HXA-GCUFF-SLB)  
M (HXA-GCUFF-MLA / HXA-GCUFF-MLB)  
L (HXA-GCUFF-LLA / HXA-GCUFF-LLB)  
XL (HXA-GCUFF-XLLA / HXA-GCUFF-XLLB)

Manufacturer  
**OMRON HEALTHCARE Co., Ltd.**  
53, Kunotsubo, Terado-cho, Muko, KYOTO, 617-0002 Japan

EU Representative  
**OMRON HEALTHCARE EUROPE B.V.**  
Scorpius 33, 2132 LR Hoofddorp,  
THE NETHERLANDS  
www.omron-healthcare.com

Importer in EU

Production facility  
**OMRON DALIAN Co., Ltd.**  
No. 3, Song Jiang Road, Economic and Technical Development Zone,  
Dalian 116600, China

**OMRON HEALTHCARE SINGAPORE PTE LTD.**  
438A Alexandra Road, #05-05/08, Alexandra Technopark,  
Singapore 119967  
www.omronhealthcare-sg.com

**OMRON HEALTHCARE, INC.**  
1925 West Field Court Lake Forest, IL 60045 U.S.A.  
OmronHealthcare.com

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
(A)


## Variety of Cuff

- Cuff tube lengths are 175 mm (7 inches) and 1,000 mm (39 inches).
- For the cuff tube length that can be connected, see the instruction manual for the OMRON non-invasive blood pressure monitors.

CUFF SIZE	Arm Circumference	
	(cm)	(inch)
SS	12 to 18	5 to 7
S	17 to 22	7 to 9
M	22 to 32	9 to 13
L	32 to 42	13 to 17
XL	42 to 50	17 to 20

## How to Wrap the Cuff

## Instructions for Use

- This product is an upper arm cuff for OMRON non-invasive blood pressure monitors.
- The cuff is not a sterilized product.
- The cuff is reusable.
- Warranty not applicable unless damaged or malfunctioning when package is first opened.
- Please report to the manufacturer and the competent authority of the Member State in which you are established about any serious incident that has occurred in relation to this device.

## Notes on Safety

### Warning

- Do not bend cuff tube during inflation and deflation, particularly after a change of body position.
- Do not wrap the cuff on the following parts:
  - An upper arm on which intravenous drip or a blood transfusion is being performed.
  - An upper arm on which an SpO<sub>2</sub> sensor, IBP catheter, or other instrument is attached.
  - An upper arm with a shunt for hemodialysis.
- If frequently performing NIBP measurement using a cuff over an extended period of time, periodically check the patient's circulation. In addition, wrap the cuff as indicated in the cautionary points.
- Ensure that the cuff tube does not become wrapped around the patient's arm or other body part.
- Only connect this cuff to the specified products. Risk of accident.

### Caution

- If a doctor has indicated that the patient has hemorrhagic diathesis or hypercoagulability, check the condition of the arm after measurement.
- Make sure the connections are tight. An air leak will prevent correct measurement.
- Use the appropriate cuff size to ensure correct measurements. If too large a cuff is used, the measured blood pressure value tends to be lower than the actual blood pressure. If too small a cuff is used, the measured blood pressure value tends to be higher.
- Do not sterilize by autoclave, ultraviolet irradiation, radiation, gas sterilization (EOG, formaldehyde gas, high concentration ozone, etc.). In addition do not clean or sterilize using a method other than that indicated in "Cleaning and Storage" in these directions, and do not wash. The cuff may become deformed and correct measurement will not be possible.
- Do not store in a location where corrosive gas is generated. The cuff may become deformed and correct measurement will not be possible.
- Do not allow any liquids inside the cuff. If a liquid gets in the cuff, dry the inside well. If used with even a slight amount of liquid in the cuff, the liquid may enter the monitor to which the cuff is connected and cause the monitor to malfunction.
- Do not inflate the cuff when it is not wrapped on the upper arm. This may damage the cuff.
- Do not leave in direct sunlight or ultraviolet light. This may damage the cuff and cause it to deteriorate.
- The cuff is a consumable product. If it becomes damaged, replace it. Correct measurement may not be possible with a damaged cuff.
- Do not use the cuff if it is damaged or has holes. It may burst during measurement.

## Description of Icons (Fig. A)

- |   |   |
|---|---|
| ① Attachment diagram  | ⑩ Need for the user to consult the instruction manual                           |
| ② Range indicator of arm circumferences to help selection of the correct cuff size. | ⑪ EU-representative   |
| ③ Range pointer and brachial artery alignment position                              | ⑫ Medical device  |
| ④ Caution   | ⑬ Need for the user to follow the instruction manual thoroughly for your safety |
| ⑤ CE Marking  | ⑭ LOT number  |
| ⑥ Manufacturer  | ⑮ Reference catalogue number  |

Product production date is integrated in the Lot number, which is placed on the product and/or sales package: the first 4 digits mean year of production, the next 2 digits mean month of production and the last 2 digits mean day of production.

## How to Use the Cuff

- Connect the cuff to the blood pressure monitor and tighten firmly. Use an OMRON specified cuff hose.
- Wrap the cuff so that the cuff tube is on the peripheral side and the "INDEX ▼ ARTERY" mark is directly over the brachial artery. Make sure the cuff does not cover the elbow joint.
- Wrap the cuff tightly so that two fingers can be inserted.
  - If the cuff is too loose, it may pinch the patient's skin when inflated and cause subcutaneous bleeding.
- Make sure "INDEX ▼ ARTERY" is within the "RANGE" indicated on the cuff surface. If not within the range, use a different cuff size.

## Cleaning and Storage

- Wipe clean on the surface of the cuff with a cloth moistened with a 70 v/v% dilution of isopropyl alcohol, or a 76.9 to 81.4 v/v% dilution of disinfection ethanol (ethyl alcohol).
- For disposal, have a company that specializes in medical waste disposal dispose of the cuff.

## Specifications

Usage environment	Temperature: 5 ~ 40°C (41 ~ 104°F)
	Humidity: 15 ~ 85% RH (no condensation)
Storage environment	Temperature: -20 ~ 60°C (-4 ~ 140°F)
	Humidity: 10 ~ 95% RH (no condensation)
Durable period (Service life)	1 year

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